

## **CellAtlas3d Installation**

Once MorphoGraphX is install, at the command line type:

```
$ MorphoGraphX --process
```

This will print the directory where MorphoGraphX plug-ins and add-ons are located. Copy the CellAtlas3D add-on to this directory:

```
$ sudo cp CellAtlas3D.so /usr/local/share/MorphoGraphX/processes
```

If you do not have sudo rights on the machine, you can put it in your home directory:

```
$ cp CellAtlas3D.so ~/.local/share/data/MorphoGraphX/MorphoGraphX/processes
```

You must create this directory (and any missing directories leading to it) if it does not already exist.